

Product Datasheet

Human CYTL1 ELISA Kit (orb779390)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	15.63-1000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q9NRR1
Sensitivity	6.5 pg/mL
Expiration Date	Please enquire.

