

Product Datasheet

Mouse IGFBP2 ELISA Kit (orb779393)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Mouse
Range	62.5-4000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P47877
Sensitivity	23.9 pg/mL
Expiration Date	Please enquire.

