

## Product Datasheet

### Human CYCS ELISA Kit (orb779413)

**Description** The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...

**Reactivity** Human

**Range** 0.32-20 ng/mL

**Tested Applications** ELISA

**Note** For research use only

**Sample Types** tissue homogenates, cell lysates, cell culture supernates and other biological fluids

**Uniprot ID** **P99999**

**Sensitivity** 0.124 ng/mL

**Expiration Date** Please enquire.

