

Product Datasheet

Rat NGAL ELISA Kit (orb779481)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Rat
Range	0.79-50 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, urine, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P30152
Sensitivity	0.33 ng/mL
Expiration Date	Please enquire.

