

Product Datasheet

Rat APOB ELISA Kit (orb779496)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Rat
Range	31.25-2000 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q7TMA5
Sensitivity	12.1 ng/mL
Expiration Date	Please enquire.

