

Product Datasheet

Rat Plg ELISA Kit (orb779526)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Rat
Range	15.63-1000 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q01177
Sensitivity	6.9 ng/mL
Expiration Date	Please enquire.

