

Product Datasheet

Animal NPY ELISA Kit (orb779577)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Bovine
Range	31.25-2000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q6RUW3
Sensitivity	9.35 pg/mL
Expiration Date	Please enquire.

