

Product Datasheet

Mouse IL34 ELISA Kit (orb779579)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Mouse
Range	7.82-500 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Uniprot ID	Q8R1R4
Sensitivity	3.3 pg/mL
Expiration Date	Please enquire.

