

Product Datasheet

Rat TGFb3 ELISA Kit (orb779695)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Rat
Range	15.63-1000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Uniprot ID	Q07258
Sensitivity	6.3 pg/mL
Expiration Date	Please enquire.

