

Product Datasheet

Animal aHSG ELISA Kit (orb779751)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Bovine
Range	6.25-400 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P12763
Sensitivity	2.45 ng/mL
Expiration Date	Please enquire.

