


Product Datasheet

AMHR2 antibody (orb77990)

Description	Rabbit polyclonal antibody to AMHR2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant fragment corresponding to 207-440aa of human AMHR2 protein
Target	AMHR2
Form/Appearance	Lyophilized powder with 1% rAlbumin and 0.02% NaN ₃ .
Storage	Shipped at ambient temperature or with ice bags if in liquid. Store immediately at -20°C upon receipt. Avoid repeating freeze-thaw cycles.
Note	For research use only
Application notes	Reconstitute in 100μL of sterile H ₂ O with 50% glycerol.
Isotype	IgG
Clonality	Polyclonal
MW	52kDa
Dilution Range	WB: 1:500-1:4,000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293



293T and MCF-7 lysate using AMHR2 antibo...