

Product Datasheet

Rat NT ELISA Kit (orb779912)

Description	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
Reactivity	Rat
Range	46.88-3000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Uniprot ID	P20068
Sensitivity	15.2 pg/mL
Expiration Date	Please enquire.

