

Product Datasheet

Human IL1F9 ELISA Kit (orb779986)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	15.63-1000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q9NZH8
Sensitivity	5.9 pg/mL
Expiration Date	Please enquire.

