

Product Datasheet

Rat EDN2 ELISA Kit (orb780117)

Description	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
Reactivity	Rat
Range	7.82-500 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Uniprot ID	P23943
Sensitivity	2.67 pg/mL
Expiration Date	Please enquire.

