

Product Datasheet

Mouse STC1 ELISA Kit (orb780120)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Mouse
Range	6.25-400 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	O55183
Sensitivity	2.63 pg/mL
Expiration Date	Please enquire.

