

Product Datasheet

Mouse GDF11 ELISA Kit (orb780373)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Mouse
Range	12.5-800 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q9Z1W4
Sensitivity	4.8 pg/mL
Expiration Date	Please enquire.

