

Product Datasheet

Rat TNFRSF1B ELISA Kit (orb780473)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Rat
Range	0.16-10 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q80WY6
Sensitivity	0.056 ng/mL
Expiration Date	Please enquire.

