

## Product Datasheet

### Rat CELSR2 ELISA Kit (orb780533)

**Description** The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...

**Reactivity** Rat

**Range** 0.79-50 ng/mL

**Tested Applications** ELISA

**Note** For research use only

**Sample Types** Tissue homogenates and other biological fluids.

**Uniprot ID** [Q9QYP2](#)

**Sensitivity** 0.24 ng/mL

**Expiration Date** Please enquire.

