

Product Datasheet

Rat AIF1 ELISA Kit (orb780548)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Rat
Range	15.63-1000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P55009
Sensitivity	6.1 pg/mL
Expiration Date	Please enquire.

