

## Product Datasheet

### Mouse ART1 ELISA Kit (orb780613)

**Description** The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...

**Reactivity** Mouse

**Range** 0.79-50 ng/mL

**Tested Applications** ELISA

**Note** For research use only

**Sample Types** Tissue homogenates and other biological fluids.

**Uniprot ID** [P52961](#)

**Sensitivity** 0.32 ng/mL

**Expiration Date** Please enquire.

