

Product Datasheet

Animal PTHrP ELISA Kit (orb780746)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Bovine
Range	15.63-1000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P58073
Sensitivity	6.3 pg/mL
Expiration Date	Please enquire.

