

## Product Datasheet

### Mouse CTSZ ELISA Kit (orb780879)

**Description** The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...

**Reactivity** Mouse

**Range** 15.63-1000 pg/mL

**Tested Applications** ELISA

**Note** For research use only

**Sample Types** Serum, plasma and other biological fluids.

**Uniprot ID** [Q9WUU7](#)

**Sensitivity** 7.1 pg/mL

**Expiration Date** Please enquire.

