

## Product Datasheet

### Human ABCC10 ELISA Kit (orb781000)

**Description** The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...

**Reactivity** Human

**Range** 0.16-10 ng/mL

**Tested Applications** ELISA

**Note** For research use only

**Sample Types** Tissue homogenates, cell lysates and other biological fluids.

**Uniprot ID** [Q5T3U5](#)

**Sensitivity** 0.061 ng/mL

**Expiration Date** Please enquire.

