



Product Datasheet

KATNAL2 antibody (orb78110)

Description	Rabbit polyclonal antibody to KATNAL2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant fragment corresponding to 263-424aa of human KATNAL2 protein
Target	KATNAL2
Form/Appearance	Lyophilized powder with 1% rAlbumin and 0.02% Na ₂ S ₂ O ₃ .
Storage	Shipped at ambient temperature or with ice bags if in liquid. Store immediately at -20°C upon receipt. Avoid repeating freeze-thaw cycles.
Note	For research use only
Application notes	Reconstitute in 100µL of sterile H ₂ O with 50% glycerol.
Isotype	IgG
Clonality	Polyclonal
MW	53kDa
Dilution Range	WB: 1:200-1:1,000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Western blot analysis of fetal brain lys...