

Product Datasheet

Guinea pig TGFb1 ELISA Kit (orb781123)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Guinea pig
Range	15.63-1000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q9Z1Y6
Sensitivity	5.8 pg/mL
Expiration Date	Please enquire.

