

## Product Datasheet

### Mouse ABCA13 ELISA Kit (orb781140)

**Description** The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...

**Reactivity** Mouse

**Range** 78.13-5000 pg/mL

**Tested Applications** ELISA

**Note** For research use only

**Sample Types** Tissue homogenates, cell lysates and other biological fluids.

**Uniprot ID** [Q5SSE9](#)

**Sensitivity** 33 pg/mL

**Expiration Date** Please enquire.

