

Product Datasheet

Mouse DBI ELISA Kit (orb781143)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Mouse
Range	0.16-10 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Uniprot ID	P31786
Sensitivity	0.065 ng/mL
Expiration Date	Please enquire.

