

Product Datasheet

Rat FSTL1 ELISA Kit (orb781152)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Rat
Range	0.16-10 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q62632
Sensitivity	0.062 ng/mL
Expiration Date	Please enquire.

