

## Product Datasheet

### Mouse PIK3C3 ELISA Kit (orb781361)

**Description** The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...

**Reactivity** Mouse

**Range** 0.79-50 ng/mL

**Tested Applications** ELISA

**Note** For research use only

**Sample Types** tissue homogenates, cell lysates and other biological fluids

**Uniprot ID** **Q6PF93**

**Sensitivity** 0.28 ng/mL

**Expiration Date** Please enquire.

