

## Product Datasheet

### Human CYP4F2 ELISA Kit (orb781425)

**Description** The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...

**Reactivity** Human

**Range** 1.57-100 ng/mL

**Tested Applications** ELISA

**Note** For research use only

**Sample Types** tissue homogenates, cell lysates and other biological fluids

**Uniprot ID** **P78329**

**Sensitivity** 0.55 ng/mL

**Expiration Date** Please enquire.

