

Product Datasheet

Human Rock2 ELISA Kit (orb781446)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	0.32-20 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	O75116
Sensitivity	0.119 ng/mL
Expiration Date	Please enquire.

