

Product Datasheet

Mouse Igk ELISA Kit (orb781596)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Mouse
Range	78.13-5000 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P01837
Sensitivity	37 ng/mL
Expiration Date	Please enquire.

