

## Product Datasheet

### Mouse DPYSL2 ELISA Kit (orb781627)

**Description** The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...

**Reactivity** Mouse

**Range** 78.13-5000 pg/mL

**Tested Applications** ELISA

**Note** For research use only

**Sample Types** tissue homogenates, cell lysates and other biological fluids

**Uniprot ID** **008553**

**Sensitivity** 27 pg/mL

**Expiration Date** Please enquire.

