

## Product Datasheet

### Mouse CD5L ELISA Kit (orb781770)

**Description** The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...

**Reactivity** Mouse

**Range** 0.32-20 ng/mL

**Tested Applications** ELISA

**Note** For research use only

**Sample Types** Serum, plasma and other biological fluids.

**Uniprot ID** [Q9QWK4](#)

**Sensitivity** 0.136 ng/mL

**Expiration Date** Please enquire.

