

Product Datasheet

Animal HSPA1A ELISA Kit (orb781839)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Bovine
Range	1.57-100 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q27975
Sensitivity	0.59 ng/mL
Expiration Date	Please enquire.

