


Product Datasheet

Human Anti-HPF4 ELISA Kit (orb781874)

Description	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
Reactivity	Human
Range	1.57-100 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates and other biological fluids.
Sensitivity	0.59 ng/mL
Expiration Date	Please enquire.



Standard curve of Human Anti-HPF4 ELISA

...

