

Product Datasheet

Mouse AGGF1 ELISA Kit (orb781890)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Mouse
Range	78.13-5000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q7TN31
Sensitivity	34 pg/mL
Expiration Date	Please enquire.

