


Product Datasheet

Rat dLK1 ELISA Kit (orb781918)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Rat
Range	0.79-50 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates and other biological fluids
Uniprot ID	O70534
Sensitivity	0.27 ng/mL
Expiration Date	Please enquire.



Standard Curve of Rat dLK1 ELISA Kit.

