

Product Datasheet

Mouse NFE2L2 ELISA Kit (orb781941)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Mouse
Range	15.63-1000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q60795
Sensitivity	5.7 pg/mL
Expiration Date	Please enquire.

