

Product Datasheet

Rat aZGP1 ELISA Kit (orb782012)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Rat
Range	15.63-1000 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q63678
Sensitivity	7.1 ng/mL
Expiration Date	Please enquire.

