

Product Datasheet

Mouse HDLBP ELISA Kit (orb782053)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Mouse
Range	0.32-20 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q8VDJ3
Sensitivity	0.117 ng/mL
Expiration Date	Please enquire.

