

Product Datasheet

Rat APLN ELISA Kit (orb782139)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Rat
Range	125-8000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q9R0R3
Sensitivity	51 pg/mL
Expiration Date	Please enquire.

