

Product Datasheet

Mouse FIL1d ELISA Kit (orb782353)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Mouse
Range	1.57-100 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q9QYY1
Sensitivity	0.57 pg/mL
Expiration Date	Please enquire.

