

## Product Datasheet

### Rat CDHE ELISA Kit (orb782409)

|                            |                                                                                                         |
|----------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| <b>Reactivity</b>          | Rat                                                                                                     |
| <b>Range</b>               | 62.5-4000 pg/mL                                                                                         |
| <b>Tested Applications</b> | ELISA                                                                                                   |
| <b>Note</b>                | For research use only                                                                                   |
| <b>Sample Types</b>        | serum, plasma, urine, tissue homogenates and other biological fluids                                    |
| <b>Uniprot ID</b>          | <b>Q9R0T4</b>                                                                                           |
| <b>Sensitivity</b>         | 25.1 pg/mL                                                                                              |
| <b>Expiration Date</b>     | Please enquire.                                                                                         |

