


Product Datasheet

General TMAO ELISA Kit (orb782472)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	All
Range	15.63-1000 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, cell lysates, cell culture supernates and other biological fluids
Sensitivity	4.93 ng/mL
Expiration Date	Please enquire.



Standard Curve
of General
TMAO(Trimethyl...

