

## Product Datasheet

### Human CDKN1A ELISA Kit (orb782591)

**Description** The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...

**Reactivity** Human

**Range** 0.16-10 ng/mL

**Tested Applications** ELISA

**Note** For research use only

**Sample Types** tissue homogenates, cell lysates and other biological fluids

**Uniprot ID** **P38936**

**Sensitivity** 0.059 ng/mL

**Expiration Date** Please enquire.

