

Product Datasheet

Animal PTX3 ELISA Kit (orb782634)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Bovine
Range	0.32-20 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, salivas and other biological fluids
Uniprot ID	Q0VCG9
Sensitivity	0.113 ng/mL
Expiration Date	Please enquire.

