

Product Datasheet

Mouse SATB2 ELISA Kit (orb782840)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Mouse
Range	31.25-2000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q8VI24
Sensitivity	13.2 pg/mL
Expiration Date	Please enquire.

