

Product Datasheet

XIAP antibody (orb783008)

Description	XIAP antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant human XIAP
Target	XIAP
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1,k
Clonality	Monoclonal
MW	55 kDa
Uniprot ID	Q60989
Dilution Range	WB=1:500-2000, IHC-P=1:20-100, IHC-F=1:20-100, ICC=1:25
Expiration Date	12 months from date of receipt.