

Product Datasheet

GAD2 antibody (orb783055)

Description	GAD2 antibody
Species/Host	Mouse
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human GAD2 (101-200/585 aa)
Target	GAD65
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	65 kDa
Uniprot ID	Q05329
Dilution Range	WB=1:200-1500
Expiration Date	12 months from date of receipt.